



Manufacturer's Information

Test No. 17-5840/1/2: DNA Database - Saliva

Each sample set was a collection of known buccal cells from three individuals (Items 1-3), provided on either GE Healthcare FTA® Indicating Micro Cards, Bode Buccal DNA Collectors™, or Puritan® Cap-Shure® swabs. Participants were requested to analyze these items using their existing protocols.

SAMPLE PREPARATION: All substrates were prepared from human buccal cheek cells suspended in sterile phosphate-buffered saline solution. The GE Healthcare FTA® Indicating Micro Cards were spotted with 120 μ l of the cell suspension, the Bode Buccal DNA Collectors™ were spotted with 400 μ l of the cell suspension, and the Puritan® Cap-Shure® swabs were spotted with 200 μ l of the cell suspension. Items 1 and 3 were created using different male donors, and Item 2 was created using a female donor. The Items were prepared at separate times and were packaged once they were thoroughly dried. Completed sample sets were stored at -20°C until shipment on January 23, 2017.

SAMPLE SET ASSEMBLY: For each sample set, all three Items (1-3) of the same substrate type were placed in a pre-labeled sample pack envelope. The sealed sample pack envelopes were then packaged in pre-labeled Heat Seal envelopes and sealed. This process was repeated until all of the sample sets were prepared.

VERIFICATION: Laboratories that conducted predistribution analysis of the samples reported consistent results.

Manufacturer's Information, continued
Test No. 17-5840/1/2: DNA Database - Saliva

Amelogenin and STR Results

Results compiled from predistribution laboratories and a consensus of at least 10 participants.

Item	D1S1656	D2S1338	D2S441	D3S1358	D4S2408	D5S818
	D6S1043	D7S820	D8S1179	D9S1122	D10S1248	D12S391
	D13S317	D16S539	D17S1301	D18S51	D19S433	D20S482
	D21S11	D22S1045	Amelogenin	CSF1PO	FGA	Penta D
	Penta E	SE33	TH01	TPOX	vWA	
1	12,16	18,19	11,14	17,17	*	11,12
	*	10,10	13,14	*	12,13	18,22
	11,12	10,13	*	13,14	14,15	*
	32.2,32.2	11,15	X,Y	10,10	20,24	11,12
	5,10	20,26.2	6,8	9,9	17,18	
2	11,14	17,18	10,11	15,15	*	11,13
	*	11,11	10,15	*	13,15	21,24
	11,12	11,11	*	12,16	13,14	*
	28,29	15,17	X,X	12,12	21,24	12,12
	7,10	15,19	9,9.3	8,10	14,16	
3	11,12	17,22	10,14	15,16	*	12,12
	*	11,12	14,15	*	13,16	19,20
	8,11	11,11	*	14,16	13,15	*
	30,31.2	16,16	X,Y	10,14	26,26	12,14
	11,16	15,29.2	6,8	8,8	14,17	

YSTR Results

Results compiled from predistribution laboratories and a consensus of at least 10 participants.

Item	DYF387S1	DYS19	DYS385	DYS389-I	DYS389-II	DYS390	DYS391	DYS392	DYS393	DYS437
	DYS438	DYS439	DYS448	DYS449	DYS456	DYS458	DYS460	DYS481	DYS505	DYS518
	DYS522	DYS533	DYS549	DYS570	DYS576	DYS612	DYS627	DYS635	DYS643	YGATAH4 Y Indel
1	*	14	11,14	13	29	24	11	13	13	15
	12	12	19	*	15	16	*	22	*	*
	*	12	13	17	18	*	*	23	10	12
3	*	16	12,12	15	30	23	10	11	13	14
	10	12	20	*	14	16	*	23	*	*
	*	12	9	19	17	*	*	21	12	11

* Results were not received by a minimum of 10 participants for the STR and YSTR loci indicated.

The information presented here is that received from the sample manufacturer. It presents details of the design specification for the test samples and/or details of how they were prepared. This information does not necessarily represent the answers that should or could be obtained from an examination of the sample. Final interpretation of the results should be deferred until the summary report is available.