



## *Manufacturer's Information*

# **Test No. 21-5051: Quantitative Drug Analysis - Methamphetamine HCl**

---

Each sample pack consisted of two items containing different concentrations of methamphetamine HCl and lactose. Participants were requested to analyze each item and report the quantitative determination of methamphetamine HCl present in the samples.

### SAMPLE PREPARATION-

The appropriate amount of methamphetamine HCl and lactose for each item were thoroughly mixed to ensure homogeneity.

ITEMS 1 and 2 (PREPARATION): For each Item, approximately 500 mg of the powder was weighed out and deposited into a glassine bag, which was folded and secured with a label. The folded glassine bag was placed into a small zip top bag and heat sealed closed. The heat sealed bag was then placed into a pre-labeled envelope.

SAMPLE PACK ASSEMBLY: One of each of the Item 1 and Item 2 envelopes was placed into a larger pre-labeled sample pack envelope.

VERIFICATION: Laboratories that conducted predistribution analysis of the samples reported consistent results that were comparable to the preparation concentrations of methamphetamine HCl. The following methods were used to examine the items: GC, GC/FID, UV.

<u>Item</u>	<u>Preparation methamphetamine HCl</u>
1	68%
2	77%

Please note that the preparation concentration is the value used for calculations during the test preparation phase and may not necessarily represent the final concentration of the samples. It is advised to wait for the Grand Mean statistics available in the Summary and Individual Reports before evaluating performance.

*The information presented here details how test samples were prepared as well as any design specifications. This information does not necessarily represent the answers that should or could be obtained from an examination of the sample(s). Final interpretation of the results should be deferred until the summary report is available.*