



Manufacturer's Information

Test No. 22-5001: Drug Analysis

Each sample pack consisted of two items. Item 1 contained 180 mg of powder consisting of 100% testosterone. Item 2 contained 300 mg of powder consisting of caffeine. Testosterone is a Schedule III controlled substance in the United States.

SAMPLE PREPARATION-

The testosterone in Item 1 was a pure powder. The caffeine in Item 2 was also a pure powder.

ITEM 1 (PREPARATION): Approximately 180 mg of the testosterone powder was weighed out and deposited into a glassine bag which was folded and secured with a label. The folded glassine bag was placed into a small zip top bag, which was heat sealed and then placed into a pre-labeled envelope.

ITEM 2 (PREPARATION): Approximately 300 mg of the caffeine powder was weighed out and deposited into a glassine bag which was folded and secured with a label. The folded glassine bag was placed into a small zip top bag, which was heat sealed and then placed into a pre-labeled envelope.

SAMPLE PACK ASSEMBLY: One of each of the Item 1 and Item 2 envelopes was placed into a larger pre-labeled sample pack envelope.

VERIFICATION: The above substances were identified by predistribution laboratories, who used the following combined list of techniques: color tests, FTIR, and GC/MS.

The information presented here details how test samples were prepared as well as any design specifications. This information does not necessarily represent the answers that should or could be obtained from an examination of the sample(s). Final interpretation of the results should be deferred until the summary report is available.