



Manufacturer's Information

Test No. 23-5804: DNA-Semen

Each sample set consisted of two known bloodstains on white fabric (Items 1 & 2), and two questioned stains on colored fabric (Items 3 & 4). Participants were requested to analyze these items using their existing protocols.

SAMPLE PREPARATION: The stains on Items 1 – 3 were prepared using human whole blood which was either drawn into citric acid preservative bags or EDTA tubes. The stain on Item 4 was a mixture prepared from human whole blood and semen. Item 1 (50 μ L) and Item 3 (50 μ L) were created using blood from the same female donor. Item 2 (50 μ L) was created using blood from a male donor. Item 4 (50 μ L) was created using a 1:1 mixture of blood from the female donor used in Item 1 and semen from the male donor used in Item 2. Stains from different sources were prepared at separate times and were packaged once they were thoroughly dried. Completed sample sets were stored at -20°C until shipment on July 03, 2023 following completion of the verification stage.

SAMPLE SET ASSEMBLY: For each sample set, all Items (1-4) were packaged into separate envelopes and then placed together in a pre-labeled sample set envelope and sealed. The sealed sample set envelopes were then packaged in pre-labeled heat seal envelopes and sealed. This process was repeated until all of the sample sets were prepared.

VERIFICATION: All predistribution laboratories confirmed the manufacturer's expected associations. Consistent allelic results were reported for all STR and YSTR loci.

Manufacturer's Information, continued
Test No. 23-5804: DNA-Semen

Amelogenin and STR Results						
<i>Results compiled from predistribution laboratories and a consensus of at least 10 participants.</i>						
Item	D1S1656	D2S1338	D2S441	D3S1358	D5S818	D6S1043
	D7S820	D8S1179	D10S1248	D12S391	D13S317	D16S539
	D18S51	D19S433	D21S11	D22S1045	Amelogenin	CSF1PO
	FGA	Penta D	Penta E	SE33	TH01	TPOX
	vWA	DYS391	DYS570	DYS576	Y Indel	
1	13,16.3	22,24	10,11	16,16	8,12	*
	10,11	11,12	11,12	18,22	11,12	11,11
	16,18	13,15	29,30	11,16	X,X	10,10
	20,23	10,11	12,15	18,26.2	7,8	6,9
	15,16	NM	NM	NM	NM	
2	13,15	17,20	10,11	17,17	11,12	*
	10,11	12,15	13,14	23,23	11,12	11,14
	12,14	13,15	31,31	15,15	X,Y	10,12
	20,22	9,13	7,12	15,24.2	7,7	8,8
	14,18	10	18	17	2	
3	13,16.3	22,24	10,11	16,16	8,12	*
	10,11	11,12	11,12	18,22	11,12	11,11
	16,18	13,15	29,30	11,16	X,X	10,10
	20,23	10,11	12,15	18,26.2	7,8	6,9
	15,16	NM	NM	NM	NM	
4-Blood	13,16.3	22,24	10,11	16,16	8,12	*
	10,11	11,12	11,12	18,22	11,12	11,11
	16,18	13,15	29,30	11,16	X,X	10,10
	20,23	10,11	12,15	18,26.2	7,8	6,9
	15,16	NM	NM	NM	NM	
4-Semen	13,15	17,20	10,11	17,17	11,12	*
	10,11	12,15	13,14	23,23	11,12	11,14
	12,14	13,15	31,31	15,15	X,Y	10,12
	20,22	9,13	7,12	15,24.2	7,7	8,8
	14,18	10	18	17	2	

YSTR Results									
<i>Results compiled from predistribution laboratories and a consensus of at least 10 participants.</i>									
Item	DYF387S1	DYS19	DYS385	DYS389-I	DYS389-II	DYS390	DYS391	DYS392	DYS393
	DYS437	DYS438	DYS439	DYS448	DYS449	DYS456	DYS458	DYS460	DYS481
	DYS518	DYS533	DYS549	DYS570	DYS576	DYS627	DYS635	DYS643	YGATAH4
2	36,38	13	16,18	13	30	23	10	11	13
	14	10	12	20	32	17	15	9	22
	39	12	12	18	17	23	23	11	12
4-Semen	36,38	13	16,18	13	30	23	10	11	13
	14	10	12	20	32	17	15	9	22
	39	12	12	18	17	23	23	11	12

* Results were not received by a minimum of 10 participants for the loci indicated.

NM - Non-Male profile, YSTR results not expected.

The information presented here is that received from the sample manufacturer. It presents details of the design specification for the test samples and/or details of how they were prepared. This information does not necessarily represent the results that should or could be obtained from an examination of the sample. Final interpretation of the results should be deferred until the summary report is available.